Digital Pocket Refractometer

PAL-Pâtissier

Cat.No. 4508



Double Scale Brix / Baumé Water Resistant **IP65**

Compact and Robust

The PAL-Pâtisser enables digital readings in both Brix and Baume scales in less than 3 seconds. Especially for syrups, the PAL-Pâtisser is an essential tool to precisely control and asses the quantity of sugar used during cooking. Plus it is trivial to clean-simply rinse the unit under a running tape.

The Brix scale has an extended range of 0.0 to 85.0% and is practical for measurements on most food products. An automatic feature to compensate for temperature effects is built into the PAL-Pâtisser, which significantly improves the accuracy of readings taken for warm/hot samples.

SAMPLE TYPES

- Syrup
- Puree
- © Fresh Fruits © Jelly
- Sorbets
- © Fruit
- © Jam

Concentrate

Technical Specifications

Scale and Range : 0.0 to 85.0% Brix

0 to 45° Baumé

Resolution : 0.1% Brix

1° Baumé

: ±0.2% Brix Accuracy

+1° Baumé

Ambient Temperature: 10 to 40°C Sample Temperature : 10 to 75°C

(Automatic Temperature Compensation [ATC])

Response Time : 3 seconds Power supply : 2AAA Protection : IP65

(Dust-tight and Protected against water jets) Dimensions : L109×W55×H31 mm, 100g

= Related Products =

PAL-1

Cat.No. 3810

Scale and Range: 0.0 to 53.0% Brix

Accuracy: ±0.2% Brix

Sample Temperature: 10 to 75°C

(Automatic Temperature Compensation [ATC]) Dimensions: L109×W55×H31 mm, 100g



The PAL-1 is ideal for measuring almost any fruit juice, beverage or food such as soups, sauces, ketchups and jams. The PAL-1 has a simple design that is robust and easy to operate.

PEN-PRO

Cat.No. 3730

Scale and Range: 0.0 to 85.0% Brix

Precision: ±0.2% Brix

Sample Temperature : 10 to 60°C

(Automatic Temperature Compensation [ATC]

Dimensions: L160×W38×H18 mm, 70g

The PEN-PRO is a unique instrument that eliminates the need for sample handling. With the scale 0.0 % to 85.0% it is capable of measuring most sample types. In addition the PEN-PRO includes a continuous measurement mode.

ATAGO CO.,LTD.

Headquarters: The Front Tower Shiba Koen, 23rd Floor TEL: 81-3-3431-1943 FAX: 81-3-3431-1 -6-3 Shiba-koen, Minato-ku, Tokyo 105-0011, Japan overseas@atago.net http://www.atago